

Brachytherapy Transfer Board


» Model:YC-DF161B

Features

- 1 ▶ The airbag can be inflated by the procedures automatically, so the operator can carry the suspended board easily;
- 2 ▶ It can be freely in the narrow space of transfer because of the unique design of the transfer board;
- 3 ▶ The retractable storage board can be used for the dirt bucket and also for the patient's belongings;
- 4 ▶ Source applicator has no any metal material;
- 5 ▶ Six Angles source applicator can help you reach the desired location quickly;
- 6 ▶ The legs supports are all carbon fiber material;
- 7 ▶ The patient's legs can be completely flat and transferred in the most comfortable manner;
- 8 ▶ Can be combine used with the vac-lock cushion and thermoplastic sheet;
- 9 ▶ The board is very solid, even after bearing 250 kg can be lifted smoothly;
- 10 ▶ Equipped with management software, it can store and record the patient's treatment information.



Parameters



Mechanical Brachytherapy Transfer Board

Model:YC-DF161MRI



Features

- 1 There are rails on the transfer bed board, which can be transferred by means of pushing, pulling;
- 2 We can provide the special material used for the MR and carbon fiber material used for the CT;
- 3 The bed has the safety lock design and the height can be adjusted;
- 4 The carbon fiber material transfer bed is operate very easily;
- 5 Seamless transport is possible between Linac, CT and the MR room.



Parameters

