

Thermoplastic Sheet



Sponges



1 Insert the sponges into the groove.





2 Press down the sponges into the groove.



3 Setup is now complete.




4 Drag out the sponges and finish the treatment.



Denture

YC-DF165



Features:

1. It is easily to cast and conform with the teeth, fix to thermoplastic mask preventing head and neck rotation, translation;
2. Providing an airway opening and tongue depressor to insure breathing during treatment;
3. Compatible with head thermoplastic mask;
4. It is disposable.

Dental putty





U-shaped Head Mask


Model: YC-M10-101AIW 2.4mm thickness, 36% perforated.



U-shaped Head Mask with Larynx Extension

Model: YC-M10-101AIWE

2.4mm thickness, 36% perforated, lengthened 55mm to larynx.



3-clamp Head Mask

Model: YC-M21-102AIW 2.4mm thickness, 36% perforated.

Model: YC-M21-102BIW 2.4mm thickness, 22% perforated.



S-shaped Head Mask

Model: YC-M10-103AIW 2.4mm thickness, 36% perforated.



T-shaped Head Mask

Model: YC-M10-103AIWT 2.4mm thickness, 36% perforated.



S-shaped Head-Neck-Shoulders Mask

Model: YC-M10-201AIW 2.4mm thickness, 36% perforated.



T-type Head-Neck-Shoulders Mask

Model: YC-M10-201AIWT 2.4mm thickness, 36% perforated.

T-type pediatric Head-Neck-Shoulders Mask

Model: YC-M10-206AIWT 2.4mm thickness, 36% perforated.



4-clamp Neck-Shoulders Mask

Model: YC-M21-202AIW 2.4mm thickness, 36% perforated.

Model: YC-M21-202BIW 2.4mm thickness, 22% perforated.



4-clamp Neck-Shoulders Mask

Model: YC-M21-203AIW 2.4mm thickness, 36% perforated.


Model: YC-M21-203BIW 2.4mm thickness, 22% perforated.



5-clamp Head-Neck-Shoulders Mask

Model: YC-M21-204AIW 2.4mm thickness, 36% perforated.

Model: YC-M21-204BIW 2.4mm thickness, 22% perforated.



4-clamp Neck-Shoulders-Breast Mask

Model: YC-M21-606AIW 2.4mm thickness, 36% perforated.

Model: YC-M21-606BIW 2.4mm thickness, 22% perforated.



4-clamp Neck-Shoulders-Breast Mask

Model: YC-M10-607AIW 2.4mm thickness, 36% perforated.




2-clamp Breast Mask

Model: YC-M21-205AIW 2.4mm thickness, 36% perforated.

Model: YC-M21-205BIW 2.4mm thickness, 22% perforated.





4-clamp Belly-Pelvis Mask

Model: YC-M21-301AIW 2.4mm thickness, 36% perforated.


Model: YC-M21-301BIW 2.4mm thickness, 22% perforated.



6-clamp Pelvis Mask

Model: YC-M21-401AIW 2.4mm thickness, 36% perforated.

Model: YC-M21-401BIW 2.4mm thickness, 22% perforated.



6-clamp Neck-Breast Mask

Model: YC-M21-609BIW 2.4mm thickness, 22% perforated.



Two-fixing Bar Belly-Pelvis Mask

Model: YC-M10-501AIW 2.4mm thickness, 36% perforated.



T-type Belly-Pelvis mask

Model: YC-M10-501AIWT 2.4mm thickness, 36% perforated.



Prone Position Pelvis Mask

Model: YC-M10-701AIW 2.4mm thickness, 36% perforated.



3-point for GKT

2.4mm thickness, 12% perforated.

Three-piece x-ray Knife Mask

Model: YC-M10-602BHW
3.2mm thickness, 3 pieces per set, 22% perforated.
can be compatible with the Brainlab Machine.



Three-piece x-ray Knife Mask

Model: YC-M10-604BIW
2.4mm thickness, 3 pieces per set, 22% perforated.



Breast Mask

Model: YC-M10-702BIW 2.4mm thickness, 22% perforated.



Limbs Mask

Model: YC-M10-608AIWU
2.4mm thickness, 36% perforated, used for arms.



Limbs Mask

Model: YC-M10-608AIWL
2.4mm thickness, 36% perforated, used for legs.



SD Type Head-Neck-Shoulders Mask

Model: YC-M10-611AIW 2.4mm thickness, 36% perforated.



SD Type Belly-Pelvis Mask

Model: YC-M10-612AIW 2.4mm thickness, 36% perforated.

